









LY-ON by Jenx

The ideal platform to build a structured yet safe and comfortable posture for children who have complex postural or medical needs



Always Supportive...



EMAIL & WEB:

sales@adaptiveimports.com www.adaptiveimports.com PHONE & FAX:

Tel: (877) 767-9462 Fax: (877) 394-3465 ADDRESS:

2744 Circleport Dr. Erlanger, KY 41018 distributors of



LY-ON by Jenx

soiling may occur.

Developed with the aid of feedback from therapists and care staff, Jiraffe are delighted to bring you the Ly-on; a real solution to the challenge of positioning children who have complex postural or medical needs.

Ly-on offers the ideal platform to build a structured yet safe and comfortable posture quickly and easily with no tools required.

The versatile base offers a wide range of height adjustment - coming low enough to enable children to engage in floor based activities and raising high enough for even the tallest of care staff to have a comfortable and safe working height when engaging with the child. 15 degrees of tilt in either direction allows for opportunities to change weight distribution to aid in pressure management. This tilt can also be utilised to provide a comfortable and safe position for a child to receive their enteral feed and can also be used for therapeutic drainage positioning.

Ly-on is designed with practicality in mind and is ideal for use with multiple children due to its hygienic wipe clean material and removable washable covers. A wide range of positioning blocks ensure that even complex body shapes can be comfortably supported - and if there is a particular shape you need, just let us know, all of our positioning blocks can be customised for an exact fit.





in and out of Ly-on safe

and dignified.

Always Configurable Accessories

As standard each base comes with:

- Four braking casters one of which is direction lockable to aid in maneuverability
- Hydraulic lift action
- Tilt-in-space mechanism
- Fold down side rails
- Four side walls in green vinyl 10" high

- End bumpers in green vinyl
- Memory foam topped lying surface (covered with grey wipe clean PU material)
- One hook fastener receptive sheet
- Ten hook tabs to hold Velcro® receptive sheet in place
- Push handle

In addition to the base we recommend using at least one positioning strap to ensure the child's safety whilst using Ly-on. Ly-on is recommended for non-ambulant children and those who have passive movement.

Additional Hook Receptive Sheet

Quilted for additional comfort and machine washable at 60 degrees. The hook fastner receptive sheet allows for precise positioning of the support blocks.



Hook Tabs

Used to secure the sheet in place, these tabs are available in packs of 10.



Side Walls

Available individually in heights of 200, 250 and 300mm - always in green vinyl.



Side Lying Strap

This strap allows for controlled positioning of the pelvis, trunk or legs. Available in three sizes (size 1, 2 & 3).



Rotation Control Strap

This strap works in a similar way to the simple strap but with the added feature of being able to be pulled to control a rotated pelvis or trunk. Requires rotation control strap anchor. Available in three sizes (size 1, 2 & 3).



Rotation control frame

Used to give a point of leverage to aid in controlling pelvic or trunk rotation. Requires accessory mount.



Accessory Mount

Simply clips onto the side bar this mount is used to receive either the rotation control anchor or activity frame. Push button release makes it easy to position this accessory anywhere on the side rails. Multiple accessory mounts can be used at any one time.



Activity Frame

Used to support either the tablet mount or tie on beam this accessory allows for the precise positioning of items for the child to engage with. Requires accessory mount. Available in two sizes (size 1 & 2).



Activity Beam

An ideal platform on which to tie toys or items to engage the child's attention whilst they are using Ly-on. Requires activity frame and accessory mount.



This accessory allows for tablets to be positioned for the child to touch or view - ideal for occupying children when using Ly-on. Requires activity frame and accessory mount. Designed to support tablets of all sizes.

Triangular Shelf Block

Comes with Velcro® loops and creates a level surface on which to place activities for the child to engage with.

Positioning Blocks

These blocks allow for the positioning and support of the child in numerous positions as required. Available in various shapes and sizes.

Bridge Block

Ideal for supporting legs in a side lying position and creating leg abduction. Available in three sizes (size 1, 2 & 3).

Roll Pillow

Allows for a neutral head position.

Flat Pillow

Come as a pair to aid in correct height setting for a neutral head position.

Accessory Box

A clever box that safely attaches to the side rails - ideal for placing enteral feed giving sets in. The accessory box is a useful storage area for numerous items that may be needed by the child.

Ly-on Accessory Cart

An ideal way to keep all of your Ly-on accessories in one safe and easy place - no more hunting around for lost cushions and straps. With four locking castors as standard the accessory cart can be easily moved around the room and is always there when you need it.





Ly-on by Jenx Technical Specifications

	Ly-on LY01-03
Max. User Height	63 in
Max. User Weight	154 lbs
Support Platform Length	66 in
Support Platform Width	22 in
Side Wall Height 10" standard (3 sizes available)	8" / 10" / 12"
Tilt-in-Space Angle	-15° to +15°
Footprint	77" x 44"
Support Platform Height from Floor (min-max)	20.5" to 40"
Weight of Product	168 lbs

Always here

To arrange a demonstration or assessment with the Jenx Ly-on, just get in touch!

My name is:

Please ring me on:

Or email me at:



Fax: (877) 394-3465

ADDRESS:

