



PLAY

ENGAGE

PARTICIPATE



Molded positioning pod incorporated into a bean bag base to maintain proper postural alignment. Ideal for use in home, school and clinical environments.

Inspired  
by drive

Pool Blue



Standard Features:

Removable Cover with Flip-up Headrest and Pelvic Belt

Positioning Pod

Bolster Support

Bean Bag Base



Starfish Red\*



Whale Blue\*



Jellyfish Purple\*



Mermaid Pink

Seahorse Green



(Shown in XS, S and L respectively)



\* Starfish Red, Whale Blue and Jellyfish Purple available in size Large and Extra Large P Pods only.



# Accessories



## Activity Frame

For use during play therapy. Allows child to participate in stimulating play while positioned in the P Pod. (Toys not included)

	Small	Medium
For Use With	XS	S / M / L / XL
Overall Width	30.25"	30.25"
Overall Length	15.75" - 34.75"	15.75" - 34.75"
Overall Height	31" - 44"	46" - 59"



## Freestanding Tray

Provides a stable surface for work and play. Height and angle adjustable.

	Small	Medium	Large
For Use With	XS	S / M	L / XL
Tray Width	20.5"	25"	29.25"
Tray Depth	21.5"	23.5"	27.25"
Overall Height	18" - 22.75"	12.25" - 32"	30" - 39.5"

## Lap Tray

Beanbag underside to stabilize during use.

	Small	Medium
For Use With	XS	S / M / L / XL
Tray Width	18.75"	20.5"
Tray Depth	12.25"	15.25"



## Mobility Base

Fits through standard door way. Cushioned handles, front locking casters and foot operated brake. Compatible with most patient lifts. (For Indoor Use Only)

	Small	Medium	Large
For Use With	XS	S / M	L / XL
Weight Capacity	88 lb	165 lb	220 lb
Overall Width	23"	26.75"	26.75"
Overall Length	40.5"	44"	48.25"
Overall Height	41.75"	42"	42"



## Foot Bolster

Adjustable secure foot support.



## Lateral Bolsters

Reduce width for increased support.



## Contoured Chest Harness

For additional upper body control.



## Safety Information

Children with challenging behavioral problems or severe involuntary movement may not be suitable for the P Pod due to such excessive body movement posing the risk of severe imbalance during use. A careful risk assessment with a healthcare professional is highly recommended.

## Measurements for P Pod



Always keep the tape measure straight across the body.  
Do not wrap the tape measure around the body.



### Hip Width

Measure from outside of the hips to the outside of the hips at the widest point.



### Trunk Width

Measure across the chest from armpit to armpit at the widest point.



### Seat Depth

Position the leg at a right angle. Measure from the buttock to the back of the knee.



### Axilla Height

Position the leg at a right angle. Measure from the rear of the buttock to the armpit.



### Shoulder Height / Head Height

Position the leg at a right angle. Measure from the bottom of the rear of the buttock to the top of the shoulder for shoulder height; to the top of the head for head height.



### Knee-to-Heel Height

Bend the hip and knee at a right angle. Measure from the back of the knee to the bottom of the heel.



MODEL NUMBER	APP-1000XS	APP-2000S	APP-3000M	APP-4000L	APP-5000XL
Product Size	Extra Small	Small	Medium	Large	Extra Large
Upholstery	MP/PB/SG	MP/PB/SG	MP/PB/SG	MP/PB/SG/WB/JP/SR	MP/PB/SG/WB/JP/SR
Hip Width	10"	10.25"	12.25"	14.25"	16.5"
Trunk Width	8.75"	10"	10.5"	12.5"	13.5"
Seat Depth	9"	11.5"	13"	13.75"	17"
Axilla Height	10.25"	11.5"	12.5"	15"	16.5"
Shoulder Height	14.5"	16.5"	17.25"	20"	21.75"
Head Height	19.75"	22"	24.5"	26"	28.75"