

Seating

Standing

Sleeping

Therapy

Standz by Jenx

A versatile standing system offering prone or supine standing that accommodates abducted or neutral leg positioning for children between 1 - 9 years.

Key features

- Tool free adjustments
- Individual leg abduction and adjustment
- Offers Prone or Supine standing
- Antimicrobial covers to aid infection control in multi-user environments

Age range: 1 to 9 years
Max user height prone/supine: 54¼" / 54¼"
Max user weight: 99 lbs
Individual leg abduction: Neutral to 30°

Always Supportive...



Standz by Jenx

A versatile standing system offering prone or supine standing that accommodates abducted or neutral leg positioning.

Inspired by the clinical benefits of abducted standing and created alongside therapists and carers all over the world ensures the perfect combination of clinical positioning with simplicity of use. Standz offers Prone or Supine standing and easily accommodates either abducted or neutral leg positioning for children between the ages of 1-9 years. Leg length discrepancy is easily accommodated by the individual control for knee and foot positioning, and all adjustments are tool free to enable quick and effective clinical postural adjustments. With an extensive range of accessories, Standz can accommodate children with a wide range of postural requirements.

INNOVATIVE FEATURE

Tool Free Adjustments

No more hunting for the right sized allen key, cleverly designed wing knobs offer a comfortable grip to enable effective clamping with minimal effort – meaning growth adjustments have never been easier.

INNOVATIVE FEATURE

Angle Gauges

Built in angle gauges for frame angle and individual leg abduction - easy to read gauges enable precise positioning and allow for the correct position to be repeatedly achieved.

INNOVATIVE FEATURE

Tray Holder (optional accessory)

This clever accessory offers convenient storage when the tray is not in use.

INNOVATIVE FEATURE

Safety Lock

Built in safety lock for frame angle adjustment – Ideal for those busy environments where there is potential for accidental activation of an unlocked mechanism.

INNOVATIVE FEATURE

Antimicrobial Covers

As standing aids the digestive system and gravity can take it's course, our covers have been designed to wipe clean, they also contain an antimicrobial agent to aid infection control.

INNOVATIVE FEATURE

Individual Leg Abduction

Standz offers the ability to independently adjust either leg from 0 – 30 degrees of abduction depending on the child's clinical needs. Simply adjusted through a wing knob clamp, this action is quick and simple to adjust for each standing session. An incorporated angle gauge makes it easy to see the degree of abduction for each leg.

INNOVATIVE FEATURE

Individual Leg Adjustment

Adjustment for both knee and footplate height – Ideal for accommodating leg length discrepancy. This tool free adjustment can be made in seconds and securely clamps into place using our clever clamping marquee design.



Standz



Standz Frame

Comes complete with:

- 4 lockable swivel casters
- Tilt mechanism with foot pedal release and safety lock
- 2 leg posts marked with numbers for easy positioning of leg supports
- 0 – 30 degree independent abduction movement of leg posts
- Individual footplates with plantar and dorsiflexion movement
- Angle gauge to read frame angle – frame angle range of 0 – 90 degrees
- Pelvic support board and cover pad in a choice of four colors
- Trunk support board and cover pad in choice of four colors

In addition to the base you will need to add the following accessories to ensure a minimum safe configuration.

Prone – sandals, knee support, pelvic support & trunk support.

Supine – sandals, knee support, knee support straps, pelvic support, trunk support & head support.

Always Configurable Accessories



Flat Pad Head Support

A simple yet effective head support pad for use in supine for children with reasonable head control. It can also be used in prone to extend the height of the chest pad.

MH09-XX

7"(w) x 2"(d) x 4"(h)



Multigrip Head Support (supine standing only)

The original and innovative Jenx Multigrip™ head support is available in 2 sizes for all Jenx seating/standing systems.

1 Link MH06 11.5"(w) x 2"(d) x 3.5"(h)

2 Link MH0101 18.7"(w) x 2"(d) x 4.7"(h)



Molded PU Head Support

For use in Prone. Simple headrest for children who require minimal head support. Comes mounted on a flexible bracket.

MH07 XX



Flat Shoulder Support

A simple yet effective wide pad for supporting the shoulders when using Standz in a supine position. Can also be used in conjunction with the pelvic support board - upper back pad in prone.

STZ 014-XX

2"(d) x 12"(w) x 3.5"(h)



Lateral Supports with Strap

Firm lateral support pads offer stability to the trunk, pelvis or both together, and the incorporated support strap can be used to facilitate extension or rotation.

STZ 012-XX

5"(d) x 1"(w) x 4"(h)



Pelvic Support Board

For use in prone positioning only, the pelvic support board offers a full contact surface to facilitate the full extension of the pelvis and prevent pelvic rotation, pelvic obliquity and unwanted hip flexion.

STZ 015-XX

2"(d) x 6"(w) x 6"(h)



Pelvic Support Board - Upper Back Pad

Used in conjunction with the pelvic support board, the back pad offers support for the child's upper back. A head support can also be mounted to the back pad to help keep a neutral head position for those children who hyper extend their neck and head when standing in prone.

STZ 016-XX

6"(w) x 2"(d) x 6"(h)



Multi Axis Knee Supports

Use for supine standing only. For use where there is the need to accommodate fixed knee flexion.

STZ 017-XX

5.5"(w) x 4.5"(h)



Standard Knee Support

Individual knee supports offering the ability to accommodate knee flexion when used in supine (requires knee straps when used in supine). Width adjustable to correctly gutter the knee in position.

STZ 007-XX

3-4"(w) x 5"(d) x 4.5"(h)



Standard Knee Strap

Intended for use in conjunction with the standard knee supports, the standard knee strap simply holds the knee in the desired position and prevents knee flexion from occurring.

MU54-XX



Sandals (Pair)

Available for Standz in either size 1 or 2, these sandals ensure solid foot placement to give the foundation for effective leg positioning and safe standing.

S16 3.9"(w) x 6.5"(d)

S2 4.3"(w) x 8.6"(d)

Accessories continued



Multi-position Tray

Cleverly designed to be used in both prone and supine, the tray is angle adjustable to ensure the correct position in either application.

STZ 008-XX
23"(w) x 15"(d)



Tray Extension Wings

Strongly recommended for use in supine and just as effective when used in prone, these wings attach to the rear of the tray to offer a support surface for the child's elbows and prevents arms.

STZ 011-XX
6"(w) x 6"(d)



Tray Holder

A clever addition that offers a solution of where to store the tray when it is not being used.

STZ 010

Standz by Jenx Technical Specification

	Units	Prone	Supine
Angle Range	degrees	0° - 90°	90° - 0°
Age Range (approx)	yrs	1 - 9 yrs	1 - 9 yrs
Max User Height	inch	54¼	54¼
Frame Width	inch	25	25
Frame Length	inch	35¼	35¼
Minimum Height <small>(footplate to top of chest pad)</small>	inch	22¾	22¾
Max Height <small>(footplate to top of chest pad)</small>	inch	39¼	39¼
Chest Pad Width <small>(min/max)</small>	inch	7-11¾	7-11¾
Hip Pad Width <small>(min/max)</small>	inch	7-11¾	7-11¾
Kneeblock Width <small>(centre to centre) (min/max)</small>	inch	4¾	4¾
Height of Stander <small>(floor to top of upper pad) (max)</small>	inch	42½	42½
Height of Stander <small>(Horizontal - floor to support board)</small>	inch	n/a	31¾
Leg Abduction	degrees	30°	30°
Max User Weight	lbs	99	99
Footplate Angle <small>(dorsiflexion)</small>	degrees	10°	-25°
Footplate Angle <small>(plantarflexion)</small>	degrees	-25°	10°



Standz Color Range



distributors of

