


Date	Buyer	PO#	Mark For
Address			
Contact	Phone ( )	Fax ( )	

**Step 1 – Seat Sizes (Required options are indicated below by † symbol)**






	ThevoTwist seat size 1			ThevoTwist seat size 2			HCPCS	Description
	£ Blue	£ Red	\$	£ Blue	£ Red	\$		
†Seat	668821	668811	\$480	678821	678811	\$520	E2294	Individual body adapted, integrated seat/back with MIS-Microstimulation®; 145° recline gas spring operated; growth adjustable seat upholstery; articulated footrest height and angle adjustable; height, depth angle adjustable swing away head rest bracket; front/rear facing seat; ventilated/washable upholstery in 2 colors
†Back	668822	668812	\$520	678822	678812	\$580	E2293	
†Recline	668823	668813	\$280	678823	678813	\$280	E1014	
†Hip sup.	668824	668814	\$200	678824	678814	\$200	E0956	
†Calf strap	6635	6635	\$35	6735	6735	\$35	E0995	

\$245	E0955	†Headrest is a required option!
\$125	E1028	†Headrest bracket is a required option!

		
6387/ 6587	6598	6599
£ †Standard	£ †Neck support	£ †Shell

Seat Size Chart	Lower Leg	Seat Depth	Seat Width	Seat to Head	Max Height	Hip Angle
Size 1	7 – 11½"	6¾ – 12"	7¾ – 11"	18 ¼ – 26 ¾"	49"	90 - 133°
Size 2	8½ – 14½"	9¼ – 15"	10 – 13½"	25 ¼ – 30 ¾"	60"	90 - 145°

**Step 2 – Frame Options (Required options are indicated below by † symbol)**

**A-chassis** - Folding aluminum frame, tilt-in-space, dual spring suspension, cargo basket, tilt-in-space lever at push bar

**Q-chassis** - hi-low indoor base/feeder/therapy frame, height/ tilt adjustable, pedal controlled gas spring, lockable rear wheels



	HCPCS	SIZE 1	SIZE 2
A-chassis	E1234	£ 6202 \$1975	£ 6402 \$2025
†Castor pin lock, pair	K0073	67830 \$92	67830 \$92
†Bus tie-down loop kit (*)	K0108	67839 \$175	67839 \$175
Pneumatic wheels, full set	K0070	£ 6646 \$205	£ 6646 \$205
Frame padding, reversible	K0108	£ 6622 \$105	£ 6622 \$105
Umbrella w/ A-chassis	K0108	£ 6594 \$10	£ 6594 \$10
Tip Assist	K0108	£ 6750 \$69	£ 6750 \$69
Cargo basket	K0108	NC	NC

	HCPCS	SIZE 1 & 2
Q-chassis	E1399	£ 6666 \$1399
Push bar for Q-chassis	K0108	£ 6660 \$275

**Q-chassis Technical Specs \*\***

Capacity	88 lbs
Seat tilt	-7 to +31°
Floor to seat	11 - 28"
Seat direction	front facing

**ThevoTwist Technical Specs \*\***

	SIZE 1	SIZE 2
Seat capacity	77 lbs	88 lbs
Frame capacity	99 lbs	110 lbs
Unit weight	36 lbs	42 lbs
Max seat tilt	31°	
Seat suspension, lockable	Yes	
Upholstery	Trevira-CS	
Foot angle	78° - 100°	
Size of foot rest	11½" x 7⅞"	13" x 8¼"
Floor to handle	30 – 46"	32 – 47"
Floor to seat	19"	19½"
Max width	23"	27"
Fold length	32"	35½"
Fold height	18½"	18½"
Brakes	foot operated	
Seat direction	front & rear facing	
Front wheels	7½" solid	
Rear wheels	10" solid	

**Z-chassis** - indoor base/feeder/therapy frame, manual height/ tilt adjustable, lockable rear wheels



	HCPCS	SIZE 1 & 2
Z-chassis	E1399	£ 6691 \$899

**Z-chassis Technical Specs \*\***

Capacity	88 lbs
Seat tilt	-15 to +35°
Floor to seat	15 - 23"
Seat direction	front facing

\*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

\*\*More details in Technical Manual

**Step 3 – Positioning & General Options (Required options are indicated below by † symbol)**

<b>Lower Body Positioning</b>		HCPCS	Unit	Small	Medium	Large	Price
Foot straps with cam buckle		E0952	Pair		£ 6648		\$30
Ankle Huggers, Bodypoint®		K0108	Pair		£ 213151		\$205
Footrest pad, height 1¾"		K0108	Each	£ 10622		£ 10023	\$95

<b>Mid Body Positioning</b>		HCPCS	Unit	Small	Medium	Large	Price
Abduction block, padded		E0957	Each	£ 6620		£ 6720	\$138
Hip supports	2 point, padded, SR buckle	E0978	Each	£ 8420		£ 8470	\$71
	4 point, padded, SR buckle	E0978	Each	£ 8421		£ 8471	\$170
	Soft pelvic harness	E0978	Each	£ 6614		£ 6714	\$139
Crotch strap (fits 8420,8470,8421,8471)		E0978	Each		£ 6341		\$25
Tray, height adjustable, clear with picture window		E0950	Each	£ 66281		£ 67281	\$190
Mounting hardware tray – †required when ordering 6628/6728		E1028	Each	£ 66282		£ 67282	\$205
Arm rest, height, depth and angle adjustable		E0973	Pair	£ 8431		£ 8478	\$375

<b>Upper Body Positioning</b>		HCPCS	Unit	Small	Medium	Large	Price
Chest harness supports	Butterfly style, neoprene	E0960	Each	£ 6385		£ 6585	\$145
	Vest style	E0960	Each	£ 6631	£ 6632	£ 6731	\$145
	H-style, posture support, AEL	E0960	Each		£ 14813		\$185
	PivotFit shoulder harness, dynamic, Bodypoint®	E0960	Each		£ 214270		\$185
	PivotFit shoulder harness, Bodypoint®	E0960	Each		£ 214250		\$165
Trunk supports	Lateral, Dynamic, MIS-Microstimulation®			£ 6383		£ 6583	\$255
	Rigid lateral	E0956	Pair	£ 6377		£ 6577	\$200
	Flexible lateral w/ chest strap	E0956	Pair		£ 6596		\$285

<b>General Accessories</b>		HCPCS	Unit	Small	Medium	Large	Price
Grab bar, padded, black		K0108	Each	£ 6621		£ 6721	\$105
Vent tray		E1029	Each		£ 110004		\$410
Oxygen tank holder, E-tank		E2208	Each		£ 665012		\$140
Straps	Secure folded seat/frame	K0108	Each		£ 667011		\$25
	Wrist "runaway" tether	K0108	Each		£ 667111		\$25
Weather protection	Canopy, sun/rain with window, black	K0108	Each		£ 663761		\$295
	Lower body windbreaker, black for 663761	K0108	Each		£ 661852		\$95
	Canopy, sun/rain with window, grey	K0108	Each		£ 637811		\$295
	Lower body windbreaker grey for 637811	K0108	Each		£ 637812		\$95
	Canopy, rain only, clear	K0108	Each		£ 6638		\$195
	Rain shield, clear for canopy 663761/ 637811	K0108	Each		£ 6739		\$109
	Rain cape, grey	K0108	Each	£ 6340		£ 6540	\$149
	Sleeping bag, summer weight, black	K0108	Each	£ 6619		£ 6719	\$295
	Sleeping bag, winter weight, black	K0108	Each	£ 6618		£ 6718	\$299
Umbrella, angle adjustable	K0108	Each		£ 6594		\$85	

\*Bus tie-down loops meet section 19 of ANSI/RESNA WC/Volume 1 – secure all wheelchairs with four SAE J2249-compliant tie-down assemblies. Positioning supports, such as pelvic belt and chest harness are NOT intended for use as occupant restraints. Use only certified occupant restraints when transporting an occupied wheelchair in a vehicle.

\*\*More details in Technical Manual